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ENGLISH • FLOWABLE LIGHT CURE COMPOSITE

Product description: It is a polymerizable material with many types of medical procedures that benefit from the restoration of cavity classes II, III, IV and V and in this way to be able to practice minimally invasive dentistry. In addition to being a material that can be used as a base to receive some restorative materials such as composite resins or dental ceramics. In the same way as its other uses, this material is beneficial for the dental organs and the dentist due to the ability to be used as a sealing material for pits and fissures thanks to its great fluidity.

Clinical benefits:

- Symptoms relief (pain and discomfort).
- Prevents future complications.
- Very High compressive strength that increase the durability and resistance of restoration.
- Rebuild and return chewing functionality
- Natural appearance and good aesthetics.
- They avoid the use of metals such as mercury
- Duration between 5 and 7 years
- It does not give off harmful components, it is biocompatible with the body.
- Good resistance to wear, fractures and tension.

Performance characteristics:

- Excellent properties for molding and sculpting.
- High Brightness.
- Available in any vita color.
- Good polishing.
- It adheres better to the teeth than other materials.
- Easy to use because it does not require the use of complicated tools for its application.
- It is suitable for all types of patients.
- Flowable Light Cure Composite can be placed on front teeth and molars
- They can be placed in blocks (5 mm increments for class II and 4 mm for class I).
- After polishing the surface to remove any extrinsic stains, determine the shades needed for the restoration before tooth preparation or rubber dam placement. A tooth that has been dried appears lighter than normal. For this reason, a shade selected on a dry tooth will be lighter than a tooth after hydration.
- For better results is recommended to use halogen or LED lamp with a wavelength between 400-500nm that activates the photoinitiators of the resin.
- For Hybrid, Microhybrid and Nanohybrid Light cure composites the exposure takes between 20-40 seconds per layer (2mm). The exposure time should be increased up to 40 seconds when the layer is thicker than 2mm or if the resin is darker.
- Due to their higher resin content, fluid composites undergo greater contraction during polymerization.
- When choosing a flowable composite, it is considered whether mechanical or aesthetic requirements take precedence. In the first case, a material with a larger filler volume is preferred, while in the second, the particle size is crucial.
- The existence of opacifiers and dyes allows improving aesthetic results with these materials.

Shades: A1, A2, A3, A3.5, B1, B2, B3, C1, C2, C3, D2, D3 Incisal and Opaque.

Intended use: To be used when a small cavity exists and does not need to increase the size of the restoration to avoid the use a hybrid composite. (Microfilling)

Indications:

- Application for minimally invasive restorations class I, II, III, IV and V, when a small cavity exists (<2mm) and does not need to increase the size of the restoration.
- Base/Liner under indirect restorations.
- As support for undermined preparations.
- Temporary repair of acrylic and/or resin elements.
- Pit and fissure sealant.

Intended user: Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

Patient population: All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must be make restorations with Medental Flowable Light Cure Composite.

Precautions: Some patients can present allergy to the components of the formula, in this case, suspend it use. Keep away from Sunlight. Keep away from rain. Keep away the reachofchildren. Keep away from heat or sparks. Temperature limit between 4°C and 24°C. Fragile. Non-Sterile. No requires refrigeration. Only for professional dental use. Room temperature before use. Not requires cold chain to transport. Do not use after date of expiry.

Contraindications: The material must not be used when the application techniques described are not possible. Composite resins are subject to irritation in soft tissues and cells in direct contact when the monomers included in their formulation are in direct contact with these tissues, causing irritation and even pulp death if they are not 100% photocured. It is contraindicated to place composite resin in any presentation on materials that contain eugenol since these inhibit polymerization. Layers greater than 3 mm should not be placed to allow correct polymerization. In the case of flow composite, the layers should not exceed 3mm to avoid fractures due to chewing compression.

Composition: Crystalline and fumed silicon dioxide, Bis-GMA, Stabilizer and Photo initiators.

Instructions for use:

1. - Restorations class I, II, III, IV, and V. Minimally invasive(<2mm):

A.TOOTH PREPARATION: Selection of the tone, correct isolation of the tooth and remove off the cavity.

B.ETCHING TOOTH SURFACE ENAMEL: Apply enamel etchant, wait 20-30 sec. Wash with water, confirm change in the look appearance like gray chalk, if not, repeat the operation until.

C.PRIMING TOOTH SURFACE: Use a soft brush to place one layer film of bonding restoration surface, remove film excess with a gentle stream air, polymerize for 10-20 sec.

D.COMPOSITE APPLICATION: After steps A, B, C are done, apply the composite from its application tip either from the syringe or from the uni-dose tip. Must be placed 1mm thickness layers and polymerized for 20-30 sec. This technique will avoid stress due to the contraction of the composite.

E.FINISHING: The final restoration needs to be polish by traditional ways; we recommended the use of silicon diamond abrasives.

2. - Base/Liner under indirect restorations:

A. TOOTH PREPARATION: Isolate the tooth, remove off the cavity, check for expose pulp and/or for a thin layer of dentin (<1mm), if it is the case, proceed to apply Calcium Hydroxide base which will protect the pulp against thermal/electrical shocks as well as antibacterial function.

B. ETCHANT: Proceed to etch the cavity for 20-30 sec. Wash with water and check the gray chalk look of the dentin/enamel, if does not appear, repeat the operation.

C. BONDING: Using a soft brush be sure to apply a film of bonding in all the surface of the cavity, thin the bond film with a gentle stream of air and polymerized for 10-20 sec.

D. APPLICATION OF THE COMPOSITE: The Flowable Light Cure Composite can be applied from the application tip of the syringe as well as the application tip of the uni-dose tip. Up to 2mm. Thickness layer of composite can be placed at one time. Polymerized for 20-30 sec.

3. - As support for undermined preparations:

A.TOOTH PREPARATION: Isolate the tooth and remove off the cavity, cement, etc.

B.ETCHANT: Apply etch being careful to cover all the internal surface of the undercut, wait 20-30 sec. and wash with water. Confirm change on the look of the etch surface into gray chalk.

C.BONDING: Apply a film of bonding over the dentin/enamel surface carefully cover all surfaces of the undercut, use a gentle stream of air to thin the bond film and polymerized for 20-30 sec.

D.COMPOSITE APPLICATION: 1mm thickness layer are recommended will avoid contraction stress.

4. - Temporary repair of acrylic and/or resin elements:

When a small fracture in an actual composite restoration in place or and excess of wearing in the case of an acrylic crown can be perform by

following the same steps as before, etchant, bonding application of the flow composite in the cracks or in the obvious wearing zone.

5. - Pit and Fissure Sealant:

A.PREPARATION OF THE TOOTH: Correctly isolate the teeth.

B.ETCHANT: Apply enamel etch and wait for 20-30 sec. wash with water and check for the gray chalk look of the enamel, if does not appear, repeat the operation.

C.BONDING: With a soft brush, apply a thin film of bonding, remove the excess with a gentle stream of air, polymerized for 10-20 sec.

D.APPLICATION OF THE COMPOSITE: Once anterior steps are done, proceed to place the composite from the application tip either the syringe or the uni-dose tip. The fluid consistency allows adequate penetration in the fissures on the occlusal face of the molars. Polymerized for 20-30 sec.

Conservation:

- The system is designed to be used at room temperature above 24°C or 75°F.
- Shelf life at room temperature is two years for unopened kits and one year for opened kit when syringe and vital caps are replaced immediately after dispensing.
- Unopened kits can be refrigerate in (4°C-40°F) to extend shelf life. If stored in refrigerator, bring to room temperature prior to use.
- Do not store materials in proximity to eugenol containing products.
- Do not expose materials to elevated temperature or intense light.

Warnings and Precautions: In case of any change on the original function of the product, stop using it and contact the manufacturer. Read carefully and follow all instructions. Before the use of these products always use rubber dam gloves and protective glasses. The etching solution contains a 37% phosphoric acid solution, which can be harmful to skin and eyes. In case of contact with the eye, wash immediately with abundant water for at least 15 minutes and if the irritation persists consult a doctor. In case of contact with skin wash the area with soap and water. The product should not be ingested or aspirated. If accidentally swallowing induce vomit immediately and get medical attention. In case of inhalation, clean the nose with water and seek medical assistance. Ask for patient to rinse his mouth immediately after treatment. The product is designed to be used at room temperature. If refrigerate, the material may be difficult to extrude.It is recommended to cover the polymerization time indicated by the manufacturer, according to the layer thickness, because a non-polymerized resin can cause irritation or allergic reactions. Requires following the dentin preparation procedure to reduce gingival recession.

Warranty: If the use of this product on a patient or anyone else causes some anomaly such as:

- It does not fulfill the intended purpose,
 - It doesn't return the functionality of the tooth,
 - It doesn't have the properties it refers to,
 - It has damage to packaging or manufacturing,
 - The instructions for use are misunderstood,
 - Any secondary reaction or anomaly.
- Please send detailed report to: Conamco, S.A. de C.V. email: atencionclientes@conamcosacv.com.mx

MEDENTAL and/or CONAMCO S.A. de C.V. warranty the product against any fault of manufacture of the material. MEDENTAL and/or CONAMCO S.A. de C.V. do not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product.

The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUADAMED) at: <https://ec.europa.eu/tools/euadamed>.

Disposal Considerations:

1. Use New and clean Nitrite Gloves.
2. The resin must be cured in a safe container for resin, either by exposing it to the sun or with a light cure lamp.
3. Once cured, it can be discarded with solid residues.

Note: Liquid or partially cured resin must be disposed of in accordance with the regulations.

Features/Claims: Dental Light Cure Composites are part of Dental Surgery and are composed of different synthetic materials that adhere to the surface of the tooth and give it a color similar to the natural one, so its use results in a practical and cosmetic restoration. The Light Cure Composites can be used for a fractured, cracked or fissured tooth, holes, cover stains, biting and chewing problems, fill cavities caused by cavities, improve the color of the teeth. The specialist must carve the surface of the tooth so that the Resin adheres to it, it is required to etch and apply adhesive later, and the product is applied in the corresponding space and molded as desired. Once the Dental Resin is placed, it hardens with the help of ultraviolet light or blue light. To take good care of the Dental Resin, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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