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# HYBRID, MICRO HYBRID AND NANO HYBRID COMPOSITE LIGHT CURE

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## ENGLISH • HYBRID, MICRO-HYBRID AND NANO-HYBRID LIGHT-CURE COMPOSITE

### Product description:

**Hybrid and Microhybrid:** Light curing resin for all types of cavity restorations, with a unique combination of barium glass fillers and Bis-GMA based resin system for anteriors or posteriors. It is an ideal filling material for anterior or posterior restorations. Its use ranges from one week to four months, but they can also be used as definitive cements, depending on the hardening method.

**Nanohybrid:**It is a composite designed to increase aesthetics and easily polishing class VIII and IV cavities in anterior teeth. It can be used in cavities where there is occlusal loading due to it high particle resistance.

### Clinical benefits:

•Symptoms relief (pain and discomfort).

•Prevents future complications.

•Very High compressive strength that increase the durability and resistance of restoration.

•Rebuild and return chewing functionality

•Natural appearance and good aesthetics.

•They avoid the use of metals such as mercury

•Duration between 5 and 7 years

•It does not give off harmful components, it is biocompatible with the body.

•Good resistance to wear, fractures and tension.

### Performance characteristics:

•Handling during placement and subsequent polishing are important aspects of the clinical

performance of these resins. Its ability to adhere to tooth structure and ease of use are key considerations.

•After polishing the surface to remove any extrinsic stains, determine the shades needed for the restoration before tooth preparation or rubber dam placement. A tooth that has been dried appears lighter than normal. For this reason, a shade selected on a dry tooth will be lighter than a tooth after hydration.

•For better results is recommended to use LED lamp with a wavelength of ≥800mW/cm².

•For Hybrid, Microhybrid and Nanohybrid Light cure composites the exposure takes between 20-40 seconds per layer (2mm).

### Variants/Shades:

Hybrid: A1, A2, A3, A3.5, B1, B2, B3, C1, C2, C3, D2, D3 Incisal and Opaques.

Microhybrid: A1, A2, A3, A3.5, B1, B2, B3, C1, C2, C3, D2, D3 Incisal and Opaques.

Nanohybrid: A1, A2, A3, A3.5, B1, B2, B3, C1, C2, C3, D2, D3 Incisal and Opaques.

Hybrid, Microhybrid and Nanohybrid Light Cure Composites Medental usually maintain their original color and do not fade easily, which is important for dental aesthetics.

### Intended use:

**Hybrid:** For aesthetics restorations of posteriors class I, II, and V.

**Microhybrid:** For restorations in anteriors and posteriors. For all cavity classes, indirect inlays, onlays and veneers.

**Nanohybrid:** For aesthetics Restorations and anteriors, Indirect Inlays, onlays and veneers, Universal Composite for aesthetics restoration.

### Indications:

**Hybrid and Microhybrid:** Anterior or posterior teeth with the presence of class I and V caries, with adequate remaining healthy dental structure.

Anterior or posterior teeth with dental fractures with no more than 50% loss of dental structure.

Anterior or posterior teeth with mild and moderate wear facets where the occlusal stress factor can be controlled.

Anterior or posterior teeth with aesthetic alterations or alterations in shape, color and position are considered to be treated with a direct surgical procedure for their correction.

### Nanohybrid:

Direct restorations on anterior teeth:

•Class III (lesions on the proximal surface).

•Class IV (injuries affecting the incisal angle).

•Class V (cervical or gum line injuries).

Aesthetic veneers:

•Nanohybrid resins are used to create dental veneers that improve the appearance of front teeth.

•Indirect inlays and onlays:

•These restorations are fabricated outside the mouth and then cemented into place.

Core reconstructions:

•To restore endodontically treated teeth (with root canals)

**Intended user:** Dental professionals only with an academic degree, experience, and training for use dental materials, it must be manipulated by a dental professional only.

**Patient population:** All patients of any age, gender, weight, or ethnicity, with permanent dentition, who, for different reasons, must be make restorations with Medental Hybrid, Microhybrid and Nanohybrid Light Cure Composite.

**Precautions:** Some patients can present allergy to the components of the formula, in this case, suspend it use. Keep away from Sunlight. Keep away from rain. Keep away the reach of children. Keep away from heat or sparks. Temperature limit between 4°C and 24°C. Fragile. Non-Sterile. No requires refrigeration. Only for professional dental use. Room temperature before use. Not requires cold chain to transport. Do not use after date of expiry.

**Contraindications:** The material must not be used when the application techniques described are not possible. Composite resins are subject to irritation or soft tissues and cells in direct contact when the monomers included in their formulation are in direct contact with these tissues, causing irritation and even pulp death if they are not 100% photocured. It is contraindicated to place composite resin in any presentation on materials containing eugenol since these inhibit polymerization. Layers greater than 3 mm should not be placed to allow correct polymerization.

**Composition:** Crystalline and fumed silicon dioxide, Bis-GMA, Stabilizer and Photo initiators.

### Instructions for use:

#### A. CLEANING:

Since cleaning and etching may cause distinct changes in the natural color of the tooth. The selection of the appropriate color should be done prior cleaning and etching tooth. Perform a prophylaxis upon surface enamel to be etched using water and pumice. You should use a prophylactic paste that contains no acid or fluorides. The soft portions of dentine and all decay should be removed. Isolate the tooth with a rubber dam or a roll of cotton. Put a protective base of Zinc Phosphate or calcium hydroxide.

#### B. CAVITY PREPARATION:

Prepare the cavity, remove the caries. If liner is required, composite filling material is compatible calcium hydroxide and Zinc Phosphate Cement. Etch the enamel surface in the apply layer of bonding prior the application of the composite.

#### C. ETCHING TOOTH SURFACE ENAMEL:

Apply etchant over the surface enamel of tooth using a cotton pellet or brush. Allow to remain for 20-30 seconds. On primary tooth high in fluoride a 90-120 seconds. If Glass Ionomer base is used, do not etch for more than 20-30 seconds. Flush area completely for 20 seconds. Do not allow patient to rinse the mouth or permit saliva to touch etched enamel surface. Airs dry completely with an oil/moisture free air source. The surface of tooth should have a dull, white frosted appearance. If not, repeat the process allowing etchant to remain on tooth for an additional 20-30 seconds.

#### D. BONDING:

Apply a thin uniform coating of Enamel Bonding Resin immediately onto etched tooth surface. Thin the resin coating on the tooth with a gentle stream of oil/moisture free air. Avoid contact with gingival or other soft tissue.

#### E. RESTORATION:

Withdraw a suitable quality of composite: apply the composite in layers of no more than 2 mm thick. To thicken and cure the composite expose to a visible light source, continue this procedure until desired thickness is obtained, ensure that each layer is cured separately, expose the total restoration to a further 60 seconds of light. High intensity light requires less time to cure the material, so you should test to find the exposure time suitable for your type of light.

#### F. FINISHING:

A smooth surface can be obtained by using strip of aluminum oxide paper and disc ranging from course to super fine. You can polish using the conventional method for acrylics.

#### Conservation:

- The system is designed to be used at room temperature above 24°C or 75°F.

- Shelf life at room temperature is two years for unopened kits and one year for opened kit when syringe and vital caps are replaced immediately after dispensing.

- Unopened kits can be refrigerate in (4°C-40°F) to extend shelf life. If stored in refrigerator, bring to room temperature prior to use.

- Do not store materials in proximity to eugenol containing products.

- Do not expose materials to elevated temperature or intense light.

**Warnings and Precautions:** In case of any change on the original function of the product, stop using it and contact the manufacturer. Read carefully and follow all instructions. Before the use of these products always use rubber dam gloves and protective glasses. The etching solution contains a 37% phosphoric acid solution, which can be harmful to skin and eyes. In case of contact with the eye, wash immediately with abundant water for at least 15 minutes and if the irritation persists consult a doctor. In case of contact with skin wash the area with soap and water. The product should not be ingested or aspirated. If accidentally swallowing induce vomit immediately and get medical attention. In case of inhalation, clean the nose with water and seek medical assistance. Ask for patient to rinse his mouth immediately after treatment. The product is designed to be used at room temperature. If refrigerate, the material may be difficult to extrude. It is recommended to cover the polymerization time indicated by the manufacturer, according to the layer thickness, because a non-polymerized resin can cause irritation or allergic reactions. Requires following the dentin preparation procedure to reduce gingival recession.

**Warranty:** If the use of this product on a patient or anyone else causes some anomaly such as:

-It does not fulfill the intended purpose,

-It doesn't return the functionality of the tooth,

-It doesn't have the properties it refers to,

-It has damage to packaging or manufacturing,

-The instructions for use are misunderstood,

-Any secondary reaction or anomaly. Please send detailed report to: Conamco, S.A. de C.V. email: [atencionclientes@conamco-sacv.com.mx](mailto:atencionclientes@conamco-sacv.com.mx)

MEDENTAL and/or CONAMCO S.A. de C.V. warranty the product against any fault of manufacture of the material. MEDENTAL and/or CONAMCO S.A. de C.V. do not take responsibility for the use of the product outside of what is specified in this instruction sheet. Being the only obligation is the replacement of defective product.

The summary of Safety and clinical Performance (SSCP) can be retrieved from the European Database on Medical Devices (EUDAMED) at: <https://ec.europa.eu/tools/eudamed>

#### Disposal Considerations:

1. Use New and clean Nitrile Gloves.

2. The resin must be cured in a safe container for resin, either by exposing it to the sun or with a light cure lamp.

3. Once cured, it can be discarded with solid residues.

Note: Liquid or partially cured resin must be disposed of in accordance with the regulations.

**Features/Claims:** Dental Light Cure Composites are part of Dental Surgery and are composed of different synthetic materials that adhere to the surface of the tooth and give it a color similar to the natural one, so its use results in a practical and cosmetic restoration. The Light Cure Composites can be used for a fractured, cracked or fissured tooth, holes, cover stains, biting and chewing problems, fill cavities caused by cavities, improve the color of the teeth. The specialist must carve the surface of the tooth so that the Resin adheres to it, it is required to etch and apply adhesive later, and the product is applied in the corresponding space and molded as desired. Once the Dental Resin is placed, it hardens with the help of ultraviolet light or blue light. To take good care of the Dental Resin, it is recommended to avoid biting on very hard things so that this material does not splinter. Also moderate the consumption of highly pigmented foods such as coffee and wine, as well as maintain good dental hygiene and use dental protector if the patient has Bruxism.

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